


NUGA

www.nuga.kr

CE



T9




User Manual

Thank you for purchasing this product.
 Before use...
 Please read the cautions for safety in the
 user manual to use this product safely.

Quality Warranty Enclosed


Tourmanium Product

 **NUGA MEDICAL CO.,LTD.**

 185, Jiraeul-ro, Jijeong-myeon, Wonju-si,
Gangwon-do, Korea

TEL : +82-33-730-0001 / FAX : +82-33-730-0007

<http://www.nuga.kr>



ENGLISH / CE H/H / B TYPE / Rev. 00
2017.04.04

 **NUGA MEDICAL CO.,LTD.**

Product Introduction

Product Components	8
Names and Functions of Product	9
Descriptions for Controller Buttons	10

How to Use Product

How to Use Product	14
How to Take Emergency Actions during Use	17

Product Storage and Care

How to Store and Clean after Use	20
Self-check Items	21
Before any A/S request	22
Check Indicator Window of Controller	23
Product Specifications	24
Product Warranty	25

Cautions for your safety

- These cautions include descriptions for the safety of users and prevention of any damage to properties among others.
- Read through the following details to use the product appropriately.
- This Manual classifies the risk level into following three degrees and the following graphical symbols remind users of cautions for matters and operations of concerns about any risk.

Interpretations on symbols



Danger

If there is high possibility of death due to non-compliance.



Warning

If there is a possibility of an injury or death due to non-compliance.



Caution

If there is a possibility of a minor injury or damage to goods due to non-compliance.

Disclaimer

- The manufacturer, importer and/or agent shall not be responsible for any problem and/or other damage to the product which has been caused from the user's fault or operation.
- Information in the user manual has been prepared on the basis of the current product specifications.
- Details in the user manual might be different from those on the product.
- The product manufacturer will be adding and supplementing new functions at any time and new technologies shall be applied to the product continuously in the future as well.
- All the standards and functions are subject to change without any notice to users to improve performance.


 **Danger**



Do not fold the product in use.
Do not fold the product when it is stored.
(Otherwise, it might cause damage to the product.)




Avoid pulling out the power cord parts or touching the power plug with the wet hands.
(Otherwise, it might cause electrical shock or fire.)



Do not use it for other usages than the described purposes of uses.
(Otherwise, it might cause damage to the product.)



Anyone using a cardiac pacemaker should consult with a doctor before any use of this product.



Do not apply washing or dry-cleaning.
(Otherwise, it might cause damage to the product or any danger of electric shocks or fire.)



Anyone other than technicians qualified for our company's A/S should not disassemble, repair or modify the product for any reason.
(Otherwise, it might cause damage to the product or electric shock and/or fire.)

 **Warning**



Do not allow children use the product alone.



Do not place any heavy object or other electric heater on this product.
(Otherwise, it might cause injuries or damage to the product.)



Do not place the product covered with a blanket during use.
(Otherwise, it might cause damage to the product or a danger of fire.)




Avoid bending the power cord excessively or pressing it down with heavy objects.



Never use other power cord with different specifications.
(Otherwise, it might cause damage to the product.)



Avoid placing inflammable substances like candle-light, cigarette light or liquors and/or vessels containing water, heavy objects or other objects on the product.
(Otherwise, it might cause fire.)




This is not a waterproof device.
Do not immerse in water or use it in the bathroom where there is moisture.

 **Caution**




Please abide by defined time and methods.




Separate the power cord while the product is not being in use for a long time.



Keep any pin, wire or other metallic pieces from stabbing or inserting to the product.
(Otherwise, it might cause damage to the product.)



Use the product in a proper range of temperature (50~60°C).
Do not use directly on skin.
(Otherwise, it might cause burns.)




Do not expose the product to the direct sunlight.
(Otherwise, it might cause damage to the product.)



Check out the applied voltage before any use and apply the power plug only to a dedicated socket.

T9



Product Introduction

Product Components	8
Names and Functions of Product	9
Descriptions for Controller Buttons	10


Product Components

i Introduction


Please check out whether all the components are unpacked before installing the product.

Please inquire from your supplier if there is any missing or damaged component.


Some specifications and/or designs are subject to change for improving the product performance.



Main assembly



Power Cord




Extension band



User manual


* Each country's power converter for will be provided.

Names and Functions of Product



Name	Function
① Power Cord	It is connected to a socket to input the power supply to the product.
② Tourmanium Heater	The heater as the main part of the product improves the heating effects by combining the Tourmanium.
③ Controller	It is used to regulate the temperature of the heating unit.
④ Tourmanium Pocket	Pocket is provided to warm your hand.

Descriptions for Controller Buttons




Name	Function
① Power LED lamp	A lamp indicating whether the device is on or off
② Temperature indicator window (FND)	Window to indicate the temperature of the product (FND) It displays the set temperature in the setting mode, and the current temperature is displayed once the desired temperature has been set. (Temperature can be set by 1°C in the 40~70°C range)
③ Power button	A button for turning the power on and off
④ Temperature control button	A button for setting the device temperature

How to Use Product

i Before using product

- Do not fold the product in use.
- Check out the voltage in use before any use.
- Anyone using a cardiac pacemaker should consult with a doctor before any use of this product.
- Use the product in a proper range of temperature.



step. 1

Connect the power cable as shown in the figure.

i Notice


T9 can be used for a user with waist size of 26-46 inches. (By using the extension band, it can be used up to the waist size of 54 inches. Please inquire if it is used for more than 54 inches.)

step. 2

By pressing the power button, the power LED lamp is turned on. (The default temperature is set at 40°C and the current temperature is displayed.)

! Caution

There is a risk of burn caused by overheating when the device is used at HIGH temperature (e.g., 70°C) for more than an hour, so please be careful.




step. 3

Press the temperature setting button to set the temperature to the desired level (40~70°C).

i Notice


The default temperature is set at 40°C. Press the temperature adjustment button (UP/DOWN) to set the desired temperature. Once the temperature reaches the set temperature, it is maintained at the set level automatically. (The temperature can be changed by 1°C. Press and hold to change the temperature quickly.)



step. 4

Set the temperature level, and place the device on the body.

How to Take Emergency Actions during Use




step. 5

Press the power button to stop the device.

i Notice

The power supply on the product is automatically turned off after 8 hours.






step. 6

Unplug the power cord from the outlet.



Caution

If a user uses the product while wearing accessories on hands, it may cause some damage to the product.




Stop using the product immediately and remove the power plug if there is any lightning.




Stop using the product immediately and contact your supplier branch or our company if you smell any burn or detect any smoke from the product.



Stop using the product immediately and consult with a doctor if any abnormality occurs to the human body.



Turn off the power switch immediately, stop using the product and contact our company if you find out any failure or power failure from the product.



Turn off the power switch, remove the power plug immediately, and then contact our company if you have poured water into the product.

i Introduction

1. Separate the power cord while the product is not in use for a long time.
2. Be aware of the following for the storage place.
 - ① Keep the product in a place with no water.
 - ② Keep the product in a place with no concern of any adverse effect such as the air pressure, temperature, humidity, ventilation, sunlight and/or air containing dust, salt, sulfur, etc.
 - ③ Be careful of safety conditions such as the slope, vibration, shock, etc. (including those during transportation).
 - ④ Do not keep the product in a place near any chemical or gas.
 - ⑤ Avoid any direct sunlight during the storage.
3. Store the product after cleaning the components.

Environmental conditions


- During operation: temperature (0~35℃), humidity (20~80%)
- Transport and storage conditions: temperature (-25~70℃), humidity (20~80%), atmospheric pressure (700~1060hPa)

How to clean



Caution

Ensure unplugging of the power plug before cleaning. Otherwise, it might cause electric shocks.



Product top

Wipe it out softly with a piece of dry cloth.



Caution

There might be risk of injuries or fire if you fail to check out the followings. Make sure that you check out the following.

Check before use

1. Is the power cord connected to an outlet? Yes No
2. Are the power LED lamp and temperature display lamp turned on when the power is turned on? Yes No
3. Is the temperature increasing well in the product? Yes No
4. Is the temperature in the product not too high? Yes No
5. Did you check out conditions for any damage to the power cord or cable conditions? Yes No

Check after use

1. Is the power supply automatically turned off after 8 hours? Yes No
2. Is the power lamp on the controller off? Yes No
3. Is the temperature display screen turned off? Yes No
4. Did you remove the power cord from the controller? Yes No
5. Is the place to keep the product exposed to the humidity dust and/or chemicals? Yes No

Before any A/S request



Caution

Ensure checking out the followings before applying for A/S if any abnormality occurs in the product.
The product might fail to operate normally even without any failure if you use it without reading the manual.

Checklist for any abnormality in the power supply

If there is no power supply on the power lamp;

- ① Check out if the power plug on the controller is connected to an outlet.
- ② Check out if the power button in the controller is pressed to turn on the product.
- ③ Check the operating voltage.

Checklist for any abnormality in the heating

If the temperature is not increasing in the product;





- ① Check out if the power button in the controller is pressed to turn on the product.
- ② Check out if the temperature is set to too weak.
- ③ Do not disassemble the product but contact your supplier or our company.

Check Indicator Window of Controller




Caution

Do not disassemble the product arbitrarily if any error occurs, but remove the power plug, turn off the power switch, and then contact your supplier.
If you disassemble the product arbitrarily, it might cause serious injuries and/or damage to the product.

Display	Descriptions	Display location
	Disconnection of temperature sensor	Temperature Indicator Window (FND)
	Short circuit of temperature sensor	Temperature Indicator Window (FND)
	Overheating in heating material	Temperature Indicator Window (FND)
	Internally overheated	Temperature Indicator Window (FND)

Product Specifications

Product name and model name	Thermal massager T9
Manufacturer	NUGA MEDICAL CO., LTD. 185, Jiraeul-ro, Jijeong-myeon, Wonju-si, Gangwon-do, Korea
Date of manufacture / number of manufacture	Marked on the product
Product weight	2.1 kg
Product size	1370 X 230 x 33 (mm)
Cautions during use	Refer to User Manual
Operating method	Refer to User Manual
Electric rating	AC 220-240 V~, 50-60 Hz
Power consumption	90 W 

Product Warranty

Company	Product name	
	Model name	
	Warranty period	2 years from purchase
	Date of purchase	
	Place of purchase	
Customer	Phone	
	Name of customer	
	Customer number	
	Address	

※ Please keep this warranty carefully. You can receive sufficient service only by presenting this warranty with the place and the date of purchase.

- NUGA Medical Co., Ltd. provides warranty of this product according to the regulations of the Consumer Protection Act.
- Contact your local sales retailer or NUGA Medical office when the product is damaged.
- The spare parts of this product are secured for 3 years after this product is discontinued.
- You must have the place of purchase written in the warranty at the time of your purchase in order to receive warranty repairs.

1. The warranty is valid until 2 years from the date of purchase and only damages caused during normal operation are subject for warranty repairs.
2. The warranty period of this product is calculated from the date of purchase, and therefore have the store of purchase record the date of purchase on your warranty when you purchase the product (When the date of purchase cannot be verified, the warranty period is added from the day after 6 months from the date of manufacture).
3. When this product is used for commercial purposes (e.g., used for sales, in abnormal environment, etc.) or for industrial purposes, the warranty period is reduced to 6 months. (including essential spare parts)
4. This warranty cannot be reissued.

※ Warranty repair is not provided in the following cases.

- When damage is caused by natural disasters (e.g., fire, earthquake, flood, and wind)
- When replacing consumptive parts that normally wear from use
- When damage is caused by an abnormal power source or faults in a connecting unit.
- Other external reasons other than faults of the product itself.
- When damage is caused by repairs or modifications done by a person who is not a certified repair personnel of NUGA Medical's local office or service center.

 **NUGA MEDICAL CO., LTD.**

<http://www.nuga.kr>

